DISPLAY

1. Weight Indicator
2. Difference in weight compared to previous reading
3. Unit of Measure: Pounds (LB) or Kilograms (KG)
4. User number
5. Information ready to be saved
6. Time and Date

FIRST TIME SETUP

Back of the Scale

1. Select pounds (LB) or kilograms (KG).
2. Remove the red tab from the battery compartment.
3. Place the scale on a flat surface. A hard surface is best for maximum accuracy.

ADJUSTING TIME AND DATE

1. Hold down the **"SET/SAVE"** button until the HOUR starts blinking.
2. Press or hold down the **"SCROLL/HISTORY"** button to choose the desired HOUR.
3. Press the **"SET/SAVE"** button. MINUTES start blinking.
4. Press or hold down the **"SCROLL/HISTORY"** button to choose the desired MINUTES.
5. Press the **"SET/SAVE"** button. "AM" or "PM" start blinking.
6. Press the **"SCROLL/HISTORY"** button to choose "AM" or "PM".
7. Press the **"SET/SAVE"** button. Now your Clock is set.

8. Press or hold down the **"SET/SAVE"** button to choose the desired MONTH.
9. Press the **"SET/SAVE"** button.
10. Press or hold down the **"SCROLL/HISTORY"** button to choose the desired DAY.
11. Press the **"SET/SAVE"** button.
12. Press or hold down the **"SCROLL/HISTORY"** button to choose the desired YEAR.
13. Press the **"SET/SAVE"** button. You have setup Time & Date.

NOTE: The scale will turn off automatically in 20 seconds if no key is pressed.

USING YOUR BODY WEIGHT MONITORING SCALE

A. MONITORING YOUR WEIGHT OVER TIME

1. Press one of the **"USER"** buttons. Display will show TIME and DATE; then it will calibrate ("*CAL*") and will show "*0.0*" with either "LB" or "KG".
2. Step on the scale and position your feet in a manner that will distribute your weight evenly over the entire scale platform. Stand still while reading the display.
3. Press the **"SET/SAVE"** button to save your reading. Display shows WEIGHT and will alternate between DATE and TIME for a moment before recalibrating ("*CAL*") and showing "*0.0*".

NOTE: The scale will turn off automatically in 20 seconds if no key is pressed.

B. WEIGHING WITHOUT SELECTING A USER BUTTON

1. Tap on the lower corner of the scale. The scale will calibrate "*CAL*", and then will display "*0.0*" and either "LB" or "KG".
2. Step on the scale and position your feet in a manner that will distribute your weight evenly over the entire scale platform. The LCD display will flash "- - -.-" several times for approximately 3 seconds before displaying your weight. Stand still while reading the display.

If you would like to save your reading, STEP OFF THE SCALE. Press the **"USER"** button where you want to save the information. Press the **"SET/SAVE"** button to save it.

NOTE: The scale will turn off automatically in 20 seconds if no key is pressed.

ERROR MESSAGES

If "*E*" appears on the upper row of the display:

- You may have stepped on the scale before it displayed "*0.0*"
- The scale cannot detect a stable weight. Stand still while weighing.
- The load has exceeded the scale's capacity (330 lb/150 kg)

If "*E*" appears on the lower row of the display:

- The difference between your previous and current reading is 100 or more (lb or kg).

CARING FOR YOUR SCALE

Periodically wipe the scale surface down with a DRY cloth to remove dust.

REVIEWING SAVED RECORDS (HISTORY)

1. Press the desired **"USER"** button.
2. Press the **"SCROLL/HISTORY"** button to scroll down the list of saved records

NOTE: The scale will turn off automatically in 20 seconds if no key is pressed.

RESETTING AN INDIVIDUAL USER'S HISTORY DATA

WARNING **This process will erase the records for the selected user.**

1. Press the **"SET/SAVE"** button and a **"USER"** button at the same time for 3 seconds. The scale will calibrate ("*CAL*"), and then will display "*0.0*" and either "LB" or "KG".

COMPLETELY RESETTING THE SCALE

WARNING **All data will be lost. Erased data cannot be recovered.**

Make sure the LCD display is turned off.

1. Open the battery compartment.
2. Using a pin or the tip of a pen, press the button inside the small hole located near the locking tab for 3 seconds. The scale's memory has been reset.

REPLACING THE BATTERIES

This unit includes two (2) CR2032 long-life 3-volt Lithium batteries. When "*Lo*" appears on the display, replace both batteries.

CAUTION **Removing the batteries will reset the Date/Time but the recorded data will not be affected.**

Battery Replacement

a) Locate the battery cover on the back of the scale. Push the locking tab and pull the battery cover up.

b) Push the small silver tab away from the battery to loosen each battery. Hold the scale at an angle to allow the battery to drop.

c) Insert two new batteries, pressing them firmly into their slot.

d) Replace the battery cover.

e) Set the Time/Date (see ADJUSTING TIME AND DATE).

IT IS NORMAL FOR YOUR WEIGHT TO VARY DURING THE DAY AND FROM ONE DAY TO ANOTHER.

CAUTION **Weight exceeding the capacity (330 lb/150 kg) may damage your scale. This product is designed for personal use only. Not for commercial or industrial use. This unit is not waterproof; avoid contact with excessive moisture.**

Your Healthy Lifestyle

Health and fitness isn't just about counting calories and cutting fat out of your diet. It is about a balanced combination of eating right and exercising: neither dieting nor exercising alone will give you the kind of results you get when they are used together. To achieve a healthy lifestyle, you must set realistic fitness goals and reduce the number of calories you take in, regardless of whether they come from foods high in protein, carbohydrates, etc.

When you go on a diet, you may lose weight, but you'll also lose beneficial muscle and bone tissue. Instead, reduce the number of calories you are taking in by eating more healthful foods, and include exercise in your regular routine to lose weight faster by burning fat. Exercise makes you feel better; it helps improve circulation, relieves stress, and speeds up your metabolism. So create a plan of exercise and diet appropriate for your needs, body type, and general health. Soon you'll see how this powerful combination makes it easier to look and feel your best.

How does extra weight affect my health?

Extra weight can put you at higher risk for some health problems such as:

- Type 2 diabetes (high blood sugar)
- High blood pressure
- Heart disease and stroke
- Some types of cancer
- Sleep apnea (when breathing stops for short periods during sleep)

Consult with your health care professional for more information.

Helpful Eating Tips

- Make sure no more than 30% of your daily calories come from fat, with no more than 10% saturated fat.
- Never skip meals! Instead, eat sensible portions at every meal and use small, low calories snacks between meals to help control hunger.
- Learn to read food labels and follow recommended serving sizes.
- Balance your calorie intake with your activity level. Eat less on low activity days than you do on days where you enjoy rigorous exercise.
- Follow the USDA's Food Pyramid as your guide to a healthy, balanced eating plan.

Source: USS, Department of Agriculture/USS, Department of Health and Human Services

Helpful Exercise Tips

- Set goals. Write down what you want to achieve in the short and long term.
- Make exercise a habit. Choose a time every day - even 30 minutes - and stick to it.
- Vary your workout. Mix up your activities to avoid burn out.
- Recruit a friend! Working out as a team keeps you both motivated.

NOTE: These materials are not intended to replace advice from your doctor or fitness professional. Please consult with your physician before beginning any fitness program or fat or weight reduction program. Individual weight loss will vary. **Health o meter® takes no responsibility for individual results, or any claim made by a third party.**

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10 YEAR LIMITED WARRANTY

Sunbeam Products, Inc. or if in Canada, Sunbeam Corporation (Canada) Limited (collectively "Sunbeam") warrants that for a period of ten years from the date of purchase, this product will be free from defects in material and workmanship. Sunbeam, at its option, will repair or replace this product or any component of the product found to be defective during the warranty period. Replacement will be made with a new or remanufactured product or component. If the product is no longer available, replacement may be made with a similar product of equal or greater value. This is your exclusive warranty. This warranty is valid for the original purchaser from the date of initial retail purchase and is not transferable. Keep the original sales receipt. Proof of purchase is required to obtain warranty performance. Sunbeam dealers, service centers, or retail stores selling Sunbeam products do not have the right to alter, modify or any way change the terms and conditions of this warranty. This warranty does not cover normal wear of parts or damage resulting from any of the following: negligent use or misuse of the warranty does not cover: Acts of God, such as fire, flood, hurricanes and tornadoes.

What are the limits on Sunbeam's Liability?

Sunbeam shall not be liable for any incidental or consequential damages caused by the breach of any express, implied or statutory warranty or condition.

Except to the extent prohibited by applicable law, any implied warranty or condition of merchantability or fitness for a particular purpose is limited in duration to the duration of the above warranty.

Sunbeam disclaims all other warranties, conditions or representations, express, implied, statutory or otherwise. Sunbeam shall not be liable for any damages of any kind resulting from the purchase or misuse of, or inability to use the product including incidental, special, consequential or similar damages or loss of profits, or for any breach of contract, fundamental or otherwise, or for any claim brought against purchaser by any other party.

Some provinces, states or jurisdictions do not allow the exclusion or limitation of incidental or consequential damages or limitations on how long an implied warranty lasts, so the above limitations or exclusion may not apply to you.

This warranty gives you specific legal rights, and you may also have other rights that vary from province to province, state to state or jurisdiction to jurisdiction.

How to Obtain Warranty Service

If you have any question regarding this warranty or would like to obtain warranty service, please call the number below, or visit our web site at www.healthometer.com and a convenient service center address will be provided to you.

In the U.S.A: 1-800-672-5625 In Canada: 1-800-667-8623

In the U.S.A, this warranty is offered by Sunbeam Products, Inc. located in Boca Raton, Florida 33431. In Canada, this warranty is offered by Sunbeam Corporation (Canada) Limited, located at 5975 Falbourne Street, Mississauga, Ontario L5R 3V8.

PLEASE DO NOT RETURN THIS PRODUCT TO ANY OF THESE ADDRESSES OR TO THE PLACE OF PURCHASE.

© 2003 Sunbeam Products, Inc. All rights reserved. Health o meter® is a trademark of Sunbeam Products, Inc. Boca Raton, FL 33431.

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Congratulations on your purchase of this Health o meter® weight monitoring scale. It is so much more than just a scale! You can monitor your weight and weight changes up to five times over a period of time. And it can monitor up to 4 different users.

The Health o meter® Weight Monitoring Scale can help you in achieving your physical best by providing body weight readings and monitoring your weight over time.

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